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## X. References

- Shome, B., Parlow, A.F. Human follicle stimulating hormone (hFSH): First proposal for the amino-acid sequence of the subunit of (hFSH) and first demonstration of its identity with the  $\alpha$ -subunit of human luteinizing hormone (hLH). *Journal of Clinical Endocrinology and Metabolism*. 39:199-202, 1974.
- Sairam, M.R., Li, C.H. Human pituitary thyrotrophin. Isolation and chemical characterization of its subunits. *Biochemical and Biophysical Research Communications*. 51:336-342, 1973.
- Carlsen, R. B., Bahl, O. P. Swaminathan, N. Human chorionic gonadotropin. Linear amino-acid sequence of the subunit. *Journal of Biological Chemistry*. 248:6810-6827, 1973.
- The Nomenclature of Peptide Hormones, IUPSC, IUB, Commission of Biochemical Nomenclature. Recommendations. *Biochimica et Biophysica Acta* 404:152-1555, 1975.
- Frachimont, P., Human gonadotropin secretion in male subjects: 439-458 in James, V. H. T., Serio M., Martini, L., (eds.). *The Endocrine Function of Human Testis*. New York, Academic Press. 1:590, 1973.
- Marshall, J. C. Investigative Procedures. *Clinics in Endocrinology and Metabolism*. 4:545-567, 1975.
- Besser, G. M., McNeilly, A. S., Anderson, D. C., Marshall, J. C., Harsoulis, P., Hall, R., Ormston, B. J., Alexandre, L., Collins, W. P. Hormonal Response to synthetic luteinizing hormone and follicle stimulating hormone-releasing hormone in man. *British Medical Journal*. 3:267-271, 1972.
- Taymor, M. L., Thompson, I. E., Berger, M. J., Patton, W. Luteinizing hormone-releasing hormone (LH-RH) as a diagnostic and research tool in gynecologic endocrinology. *American Journal of Obstetrics and Gynecology*. 120:721-732, 1974.
- Franchimont, P., Becker, H., Valcke, J. C., Scellens, C. M., Demoulin, A., Thys, O., Bourguignon, J. P., Legros, J. J. Action of luteinizing hormone-releasing hormone (LH-RH) on gonadotrophins in normal and Pathological conditions. In Vokaer, R., DeBock, G., (eds.), *Reproductive Endocrinology*, Oxford, Pergamon. 185:202-252, 1975.
- Shaw, R. W., Butt, W. R. London, D. R., Marshall, J. C. The estrogen provocation test: a method of assessing the hypothalamic-pituitary axis in patients with amenorrhoea. *Clinical Endocrinology (Oxford)*. 4:267-276, 1975.
- Krieger, D. T., Ossowski, R., Fogel, M., Allen, W. Lack of circadian periodicity of human serum FSH and LH levels, *Journal of Clinical Endocrinology metabolism*. 34:619-623, 1972.
- John R. Marshall. Introduction of ovulation Clin. Obstet. Gyn. 21:147, 1978.

Rev: Jul.2020