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- XI. Literatur**
- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic.* Arch Inter Med. 144:1966-68, 1984.
 - Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range.* Chest 102(Suppl): 312S-316S, 1992.
 - NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline.* NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
 - Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time.* Am J Clin Path. 81: 618-622, 1984.
 - Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effect of Drugs on Clinical Laboratory Tests.* Clin Chem 18:1041, 1972.
 - Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy.* Arch Inter Med. 146:462-472, 1986.
 - Palaereti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise. Thromb Haemostasis* 58:905-910, 1987.
 11. Daten in der Datei 510(K).

	Bestellinformationen	
Bestell- Nr.	Bezeichnung	Inhalt
100356	Thromboplastin-D	10 x 4 mL
100357	Thromboplastin-D	10 x 10 ml
100352	Thromboplastin-D	10 x 2 mL
100244	Rekonstitutionsflüssigkeit	10 x 10 ml

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ENGLISH

- elements.
- All labware must be clean and free of trace amounts of detergents.
- Always follow instrument manufacturer’s instructions for proper maintenance.

Interfering Substances

- Sodium oxalate, EDTA, and heparin are not suitable anticoagulants.
- The PT may be prolonged by substances such as oral contraceptives, corticosteroids, EDTA, asparaginase, clofibrate, erythromycin, ethanol, tetracycline, and anticoagulants such as heparin and warfarin.⁵
- The PT may be shortened by substances including antihistamines, butabarbital, caffeine, oral cotraceptives, phenobarbital and vitamin K.⁵

IX. Expected Values

In multi-center studies Thromboplastin-D was used to assay a normal patient population and yielded the following results:⁶

Instrument	PT Mean(secs)	Range (+/-2SD)	N
Stago STA	11,3	10,4-12,2	129
Amelung KC 4A™	11,8	9,8-13,7	30
MLA™ Electra 1000C™	12,2	10,4-14,1	30

NOTE: The STA normal range was established using fresh samples, all other instruments with frozen samples. These values should be used as a guideline. Each laboratory should establish a Normal Reference Range (NRR) using instrumentation, blood collection methods, and testing techniques used in that laboratory. The NRR should be reestablished or at least verified when changing lot numbers of the same reagent.^{3,6} A new NRR should be established with any change in reagents, instrumentation, blood collection techniques, or anticoagulant.

The clotting time of abnormal plasmas will depend on the ISI of the reagent lot in use.

X. Performance Characteristics

Precision: Precision of Prothrombin Time results is dependent on many factors, such as the instrument, technique and the reagent in use. Thromboplastin-D precision was assessed by testing a normal and abnormal plasma on several different instruments. The results are summarized in the following table.⁹

Summary of Within-run Precision Studies, %CV (N = 20)

Sample	BioMerieux MDA®	MLA Electra 1000C	Amelung KC 4A
Normal	0,8%	0,6%	2,0%
Abnormal	0,9%	1,8%	2,1%

Sensivity: Thromboplastin-D detects deficiencies in the extrinsic pathway as determined by the Prothrombin Time test. Factor sensitivity testing was performed by diluting pooled normal plasma with factor deficient plasmas such that the final factor concentration ranged from 0-100%. PT testing of these samples was performed on the MLA-1000C instrument.¹⁰

Factor Sensitivity Testing

	Prothrombin Time (secs)			
% Factor	Factor II	Factor V	Factor VII	Factor X
100	10.9	11.3	11.2	11.4
50	10.5	12.6	12.2	13.0
40	10.5	13.1	12.8	13.3
30	10.9	13.8	13.5	14.5
20	11.4	15.2	14.3	15.9
10	13.4	17.6	16.0	19.6
0	34.0	55.9	25.9	108.8

Correlation: Correlation studies were performed against a competitor rabbit brain thromboplastin reagent on the MLA-1000C instrument.¹¹

Thromboplastin D vs. Reagent A, N = 101	PT Correlation <p>R = 0.97 y = 0.7406x + 2.945</p>	INR Correlation <p>R = 0.96 y = 0.9795x - 0.0029</p>
---	--	--

XI. References

- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic.* Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range.* Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline.* NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time.* Am J Clin Path. 81: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effect of Drugs on Clinical Laboratory Tests.* Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy.* Arch Inter Med. 146:462-472, 1986.
- Palaereti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise. Thromb Haemostasis* 58:905-910, 1987.
11. Data found in 510 (K) file.

	Ordering Information	
Cat. No.	Description	Contents
100356	Thromboplastin-D	10 x 4 mL
100357	Thromboplastin-D	10 x 10 mL
100352	Thromboplastin-D	10 x 2 mL
100244	Reconstitution Fluid	10 x 10 mL

Fisher Diagnostics® Limited Warranty

Fisher Diagnostics (FD) warrants to the purchaser only that FD products will perform as described on their labeling and product literature. Purchaser must determine the suitability of FD products for their specific applications. FD’s sole obligation will be, at its option, to either replace a non-conforming or defective product, or return the purchase price. FD DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE. Neither FD nor its affiliates shall, in any event, be liable for incidental or consequential loss or damage.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Amelung KC 4A is a registered trademark of Trinity Biotech. MLA Electra 1000C is a registered trademark of Instrumentation Laboratory.

Pacific Hemostasis® Thromboplastin-D

I. Intended Use

Pacific Hemostasis Thromboplastin-D is intended for use in performing the one-stage prothrombin time (PT) test and PT-based factor assays.

II. Summary and Principles

The PT is used as a screening tool and as a quantitative test for coagulation factors in the extrinsic and common pathways. This test will be prolonged in patients with acquired or congenital disorders that reduce the activity of factors I (fibrinogen), II (Prothrombin), V, VII, and X. The PT is also widely used to monitor oral anticoagulant therapy.^{1,2} Oral anticoagulants reduce the activity of vitamin-K dependent clotting factors (II, VII, IX, X, Protein C, and Protein S), and the PT is prolonged as a result.

The one-stage PT measures the clotting time of plasma after adding a source of tissue factor (thromboplastin) and calcium. The recalcification of plasma in the presence of tissue factor generates activated Factor Xa (F Xa). F Xa in turn activates Prothrombin to thrombin, which converts fibrinogen to an insoluble fibrin clot.

III. Reagent

Für *in vitro* diagnostic use.

Composition: <0,9 % rabbit brain tissue, 0,08 % sodium azide, 2% buffers, salts and stabilizers

Store unopened vials at 2-8°C. Reconstitute with preservative-free distilled/deionized water according to the Thromboplastin-D vial label, swirl gently and let the vial stand undisturbed for 15 minutes at room temperature. Do not invert the vial or mix vigorously. Reconstitution fluid is available if the quality of water is questionable. After reconstitution, the reagent when stored stoppered is stable for 7 days at 2-8°C, 8 hours at 37°C. Store at 2-8°C when not in use. **Do not freeze.**⁴

Mix gently before each use. Provide some mechanism, such as a magnetic stirrer, to maintain adequate suspension during use. Lack of vacuum in vials, erratic results, quality control values outside established ranges, or product color variations could indicate deterioration. However, poor performance could also be due to other factors within the test system.

Warning: Thromboplastin-D contains sodium azide. Sodium azide under acidic conditions yields hydrazoic acid, an extremely toxic compound. Azide compounds should be diluted with running water before being discarded. Upon disposal, azide compounds should be flushed with large volumes of water. These precautions are recommended to avoid deposits in metal pipes in which explosive conditions may develop.

CAUTION: This product contains animal source material. Handle and dispose of this product as if it were potentially infectious.

IV. Specimen Collection

3.2% (0.109M) trisodium citrate anticoagulant is recommended for coagulation assays. Avoid hemolysis and contamination by tissue fluids. Samples that have less than 90% of the expected fill volume should be rejected. Centrifuge blood for 15 minutes at 1500 x g. Test within 2 hours if samples are held at 22-24°C. If testing is not completed within 24 hours, plasma should be frozen at –20°C for up to two weeks or –70°C for up to 6 months. For more details on specimen collection and storage, see NCCLS Document H21-A4.³

- Do not delay mixing the blood with anticoagulant.
- Avoid foaming the specimen.
- Use only plastic or siliconized glass containers.
- Turbid, icteric, lipemic, or hemolyzed specimens may generate erroneous results.
- Freezing and thawing plasma that contains residual cells will generate damaged cell membranes that can affect results.
- Acute inflammatory reactions can shorten PT results because of elevated fibrinogen.
- Plasma samples with hematocrits outside the range of 20-55% may be improperly anticoagulated and should be adjusted appropriately.

V. Test Procedure

Materials Provided: Thromboplastin-D Reagent, 10 x 4 mL, 10 x 10 mL, or 10 x 2 mL

Materials Required, But Not Provided:

Distilled or Deionized water, or Pacific Hemostasis Reconstitution Fluid
Stopwatch or timer

Precision pipette: 0.1 and 0.2 mL

Normal and abnormal controls such as Pacific Hemostasis Coagulation Control Plasmas, Level 1, 2, and 3

Thromboplastin-D is suitable for use with manual, mechanical, photo-optical, nephelometric, or other means of clot detection. Follow manufacturer’s recommendations for proper use

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XI. Referencias

- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic.* Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range.* Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline.* NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time.* Am J Clin Path. 81: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effect of Drugs on Clinical Laboratory Tests.* Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy.* Arch Inter Med. 146:462-472, 1986.
- Palaeretti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise. Thromb Haemostasis* 58:905-910, 1987.
11. Datos encontrados en el fichero 510 (K).

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XI.

- Références**
- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic.* Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range.* Chest 102(Suppl): 312S-316S, 1992. NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline.* NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time.* Am J Clin Path. 81: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effects of Drugs on Clinical Laboratory Tests.* Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy.* Arch Inter Med. 146:462-472, 1986.
- Palaeretti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise. Thromb Haemostasis* 58:905-910, 1987.
11. Données se trouvant dans le dossier 510(K).

	Pour commander	
Réf. Cat.	Description	Contenu
100356	Thrombolastine D	10 x 4 ml
100357	Thrombolastine D	10 x 10 ml
100352	Thrombolastine D	10 x 2 ml
100244	Liquide de reconstitution	10 x 10 ml

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- XI. Βιβλιογραφία**
- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic.* Arch Inter Med. 144:1966-68, 1984.
 - Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range.* Chest 102(Suppl): 312S-316S, 1992.
 - NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline.* NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
 - Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time.* Am J Clin Path. 81: 618-622, 1984.
 - Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effect of Drugs on Clinical Laboratory Tests.* Clin Chem 18:1041, 1972.
 - Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy.* Arch Inter Med. 146:462-472, 1986.
 - Palareeti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results.* *The Bologna Oral Anticoagulant Control Exercise.* *Thromb Haemostasis* 58:905-910, 1987.
 - Τα δεδομένα βρίσκονται στο αρχείο 510(K).

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
- XI. Bibliografia**
- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic.* Arch Inter Med. 144:1966-68, 1984.
 - Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range.* Chest 102(Suppl): 312S-316S, 1992.
 - NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline.* Documento NCCLS H21-A4. NCCLS, Wayne, PA, 2003.
 - Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time.* Am J Clin Path. 81: 618-622, 1984.
 - Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effects of Drugs on Clinical Laboratory Tests.* Clin Chem 18:1041, 1972.
 - Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy.* Arch Inter Med. 146:462-472, 1986.
 - Palareeti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results.* *The Bologna Oral Anticoagulant Control Exercise.* *Thromb Haemostasis* 58:905-910, 1987.
 - Dati ottenuti dal file 510 (K).

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EC REP
 WMDE
 Bergerweg 18
 6085 AT Horn
 The Netherlands

Fisher Diagnostics
 a division of Fisher Scientific Company, LLC
 a part of Thermo Fisher Scientific Inc.
 Middletown, VA 22645-1905 USA
 Phone: (800) 528-0494
 (540) 869-3200
 Fax: (540) 869-8132

JL840799 (R2)

Thermo
 SCIENTIFIC