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- Yan Y, Chang L, Wang L., Laboratory testing of SARS-CoV, MERS-CoV, and SARS-CoV-2 (2019-nCoV): Current status, challenges, and countermeasures. Rev Med Virol. 2020 Apr 17:e2106. doi: 10.1002/rmv.2106
- Letko M, Marzi A, Munster V. Functional assessment of cell entry and receptor usage for SARS-CoV-2 and other lineage B betacoronaviruses. Nat Microbiol. 2020;5(4):562-569
- Wrapp D, Wang N, Corbett KS, et al. Cryo-EM structure of the 2019-nCoV spike in the prefusion conformation. Science 2020;13:367 (6483):1260-1263
- Hoffman M, Kleine-Weber H, Schroeder S, et al. SARS-CoV-2 Cell Entry Depends on ACE2 and TMPRSS2 and Is Blocked by a Clinically Proven Protease Inhibitor. Cell 2020;16:181 (2): 271-280.e8
- Long Q-X, Liu B-Z, Deng H-J, et al. Antibody responses to SARS-CoV-2 in patients with COVID-19. Nat Med 2020 Apr 29
- Tan C-W, Chia W-N, Qin X, Liu P, Chen M-I-C, Tiu C. A SARS-CoV-2 surrogate virus neutralization test based on antibody-mediated blockage of ACE2-spike protein-protein interaction. Nat Biotechnol. 2020;38:1073-1078
- Lou B, Li T-D, Zheng S-F, et al. Serology characteristics of SARS-CoV-2 infection since exposure and post symptom onset. Eur Respir J 2020 May 19;2000763
- Zhao J, Yuan Q, Wang H, et al. Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019. Clin Infect Dis 2020 Mar 28
- Tuillon E, Bollore K, Pisoni A, et al. Detection of SARS-CoV-2 antibodies using commercial assays and seroconversion patterns in hospitalized patients. Journal of Infection 2020 Jun 3
- US Food and Drug Administration, 2020, Accelerated Emergency Use Authorization for SARS-CoV-2 RBD IgG Test. Device Description and Test Principle <https://www.fda.gov/media/139053/download>
- Shibo Jiang, Christopher Hillyer, Lanying Du. Neutralizing Antibodies against SARS-CoV-2 and Other Human Coronaviruses. Trends in Immunology, 2020 <https://doi.org/10.1016/j.it.2020.03.007>
- National Committee for Clinical Laboratory Standards. Evaluation of Stability of *In Vitro* Diagnostic Reagents; Approved Guideline – First Edition. NCCLS Document EP25-A. 2009.
- National Committee for Clinical Laboratory Standards. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition. NCCLS Document EP5-A2. 2004.

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