



Azma Elixir Pajoo

Azma cDNA Synthesis Kit

DESCRIPTION

The Azma cDNA Synthesis Kit is a complete system for efficient synthesis of first strand cDNA from mRNA or total RNA templates. The kit uses Thermo-resistant H Minus M-MuLV Reverse Transcriptase, which has no RNase H activity. Therefore, degradation of RNA does not occur during first strand cDNA synthesis, resulting in higher yields of full-length cDNA from long templates compare to other reverse transcriptases. The enzyme maintains active at 42-50°C and is suitable for synthesis cDNA from RNAs having a high degree of secondary structure.

This cDNA synthesis kit is readily compatible with various cDNA-dependent downstream applications.

Components of AZMA cDNA Synthesis Kit	Cat. No. AECD1250-50
RT Enzyme	25 μ l
5X RT Buffer	125 μ l
Oligo dT Primer (ready to use)	50 μ l
Random Hexamer Primer (ready to use)	50 μ l
dNTP Mixture (ready to use)	50 μ l
Nuclease Free Water	200 μ l

STORAGE

All components of the kit should be stored at **-20°C**.

IMPORTANT NOTES

- **Avoiding ribonuclease contamination**

RNA purity and integrity are essential for synthesis of full-length cDNA. RNA can be degraded by RNase A, which is a highly stable contaminant found in any laboratory environment.

General recommendations to avoid RNase contamination:

- DEPC-treat all tubes and pipette tips to be used in cDNA synthesis or use nuclease-free labware.
- Wear gloves when handling RNA and all reagents, as skin is a common source of RNases.
- Keep all kit components tightly sealed when not in use. Keep all tubes tightly closed during the reverse transcription reaction.
- **Make Master mix A and B in a work station and in a separate area from Extraction Room.**
- **No template negative control (NTC) is important to assess for reagent contamination.**
The NTC reaction contains every reagent for the reverse transcription reaction except for RNA template.
- **PRODUCT USE LIMITATION**
This product is developed, designed and sold exclusively for research purposes and *in vitro* use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.



PROTOCOL

Note: After thawing, gently vortex and briefly centrifuge the components of the kit. Store on ice.

- 1) Make **Master Mix A** with the following order in a **2 ml** DNase and RNase free tube:

A1 Add **1 μ l** of **Oligo dT primer** per sample.
A2 Add **1 μ l** of **Random Hexamer primer** per sample.
A3 Add **1 μ l** of **dNTP Mixture** per sample.

- 2) After mix and briefly centrifuge the **Master Mix A**, aliquot **3 μ l** of it per sample into **0.5 ml** DNase and RNase free **reaction tubes**. Keep on ice until used.
3) Make **Master Mix B** in a **new 2 ml** DNase and RNase free tube with the following order:

B1 Add **0.5 μ l** of **RT enzyme** per sample.
B2 Add **2.5 μ l** of **5X RT Buffer** per samples.
B3 Add **4 μ l** of **Nuclease Free Water** per sample.

Note: Keep **Master Mix B** on ice until used.

- 4) **Add 10 μ l** of **total RNA** or **mRNA** to every **0.5 ml reaction tube** (contained Master Mix A).
5) **Incubate** at **65°C** for **5 minutes** and then incubate at **4 °C** for **1 minute**. Briefly spin down the mixture.
6) **Add 7 μ l** of **Master mix B** to every **0.5 ml reaction tube**.
7) **Incubate** at **47 °C** for **60 minutes** and then at **80 °C** for **8 minutes**.
8) Briefly spin down the mixture and store at **-20 °C**.

Note: The synthesized cDNA can be directly used for downstream application or stored at **-20 °C**.

Note: The volume of total cDNA would be **20 μ l**.

