



## PCR AZMA Master Mix

Azma Elixir Pajoo

**Description:** PCR AZMA Master Mix is a premixed, ready-to-use solution containing Taq DNA polymerase, dNTPs, MgCl<sub>2</sub> and reaction buffers at optimal concentrations for efficient amplification of DNA templates by PCR.

**Storage Conditions:** Product may be stored at **-20°C** for up to **8 months**. Minimize the number of freeze and thaw cycles by storing in working aliquots. Mix well prior to use.

### Protocol

1. Thaw the **PCR AZMA Master Mix** at room temperature. Vortex the Master Mix and then spin it briefly in a microcentrifuge to collect the material in the bottom of the tube.

2. Prepare the following reaction mixes on ice (For a **20 µl** reaction volume)

Component	Volume	Final Conc.
2X PCR Master Mix	10 µl	1X
Forward Primer (10 pmol/µL)	1 µl	0.5 pmoles/µL
Reverse Primer (10 pmol/µL)	1 µl	0.5 pmoles/µL
DNA template	Variable	10 fg~1 µg
Nuclease-Free Water	to 20 µl	-

### General Guidelines for Amplification by PCR for under 3 Kb target sequence

#### A. Denaturation

- Generally, a **5-minute** initial denaturation step at **94°C** is sufficient.
- Subsequent denaturation steps will be between **30 seconds** and **1 minute**.

#### B. Annealing

- The annealing step is typically **30 seconds** to **1 minute**.

#### C. Extension

- Allow approximately **1 minute** for every **1kb** of DNA to be amplified.
- A final extension of **5 minutes** at **72°C** is recommended.

#### D. Refrigeration

- If the thermal cycler has a refrigeration cycle, the cycling reaction can be programmed to end by holding the tubes at **4°C** for several hours.

#### E. Cycle Number

- Generally, **25–30** cycles result in optimal amplification of desired products.
- Occasionally, up to **40** cycles may be performed, especially for detection of low-copy targets.

