



Azma Elixir Pajoooh

# RNrich Blood Kit

## Introduction

The **RNrich** Blood Kit provides all of the reagents necessary to extract RNA from whole blood. RNA purified with this kit is suitable for a variety of applications, including amplification and digestion with restriction endonucleases.

### RNrich Blood Kit components\* Cat. No. AEBRX-1104-050

<b>TR Buffer</b>	<b>20 ml</b>
<b>Activator Reagent**</b>	<b>1 ml</b>
<b>VR Buffer</b>	<b>15 ml</b>
<b>PE Buffer</b>	<b>15 ml</b>
<b>Wash Buffer(conc.) **</b>	<b>15 ml</b>
<b>Elution Buffer</b>	<b>3 ml</b>
<b>Manual</b>	<b>1</b>

\* 1.5 ml microcentrifuge tube, **molecular biology** grade water, and **absolute ethanol** are needed but are not included.

\*\* Please refer to reminder, **Activator Reagent Preparation** and **Wash Buffer Preparation** before using this kit.

### Chemical Hazard

Always wear gloves and practice standard safety precautions while using the kit. Do NOT disinfect extraction waste in solutions containing **bleach** or any other form of acid. To clean any items contaminated with the reagent, simply soak in detergent and water to remove all traces of contamination before cleaning with bleach or acidic solutions.

### Reminder

Pre-set heater block at **65°C**.

Prepare **Activator Reagent** immediately prior to use. Prepared **Activator Reagent** must be kept at **4°C**.

During RNA extraction, **never** open the microtube cap outside the laminar hood.

### Activator Reagent Preparation

Add **1 ml** of molecular biology grade water to **Activator Reagent**, and vortex it well.

**Note:** Mark the check box on the bottle.

**Note:** For the best results, the prepared Activator Reagent should be used immediately. Prepared Activator Reagent can be stored for up to 3 months or 6 months at 4°C and -20°C, respectively.

### Wash Buffer Preparation

Add **35 ml** **molecular grade absolute ethanol** to each **Wash Buffer** bottle before first use and mark the check box on the bottle.

## PROTOCOL

### Step 1: Sample Preparation

- Transfer **100 µl** of fresh whole blood sample to a 1.5 ml microcentrifuge tube.  
**Note:** The blood sample should be homogenized completely before use.

### Step 2: Tissue Digestion

**Note:** Refer to reminder before use activator reagent.

- Add **400 µl** of **TR Buffer** and **20 µl** of prepared **activator reagent** to the sample tube and vortex.
- Incubate at **65°C** for **20** minutes.  
**Note:** During incubation time, vortex vigorously the sample tube every 10 minutes.
- Incubate at **85°C** for **10** minutes.

### Step 3: Lysis

- Add **300 µl** of **VR Buffer** to the sample tube and mix by vortex.
- Keep at **room temperature** for **10** minutes and **invert** each 5 minutes.
- Centrifuge at **4°C** and **11500 g** for **10** minutes.
- Carefully transfer **300 µl** of **supernatant** to a new **1.5** microcentrifuge tube and do not disturb the phases.

### Step 4: RNA Precipitation

- Add **300 µl** of **PE Buffer** to the sample tube and **invert** slowly for 15 times.
- Incubate the samples at **-20°C** for **60** minutes.  
**Note:** In sample with low amount of RNA, incubation time can be up to overnight.
- Centrifuge at **4°C** and **12000 g** for **20** minutes.
- **Decant** the supernatant into waste.  
**Note:** While decanting, be careful not to lose the RNA pellet.

### Step 5: Wash

**Note:** Refer to reminder before use **Wash Buffer**.

- Add **500 µl** of **Prepared Wash Buffer** to the sample tube, invert slowly for 10 times and then centrifuge at **4°C** and **12000 g** for 10 minutes.
- **Decant** the supernatant into waste.
- **Repeat Step 5.**
- **Decant** the supernatant into waste, and let stand until dry.

### Step 6: RNA Elution

- Add **50 µl** of **Elution Buffer** to the sample tube.
- Resuspend the RNA for **3** minutes at **room temperature**, and then centrifuge at **10000 g** for **1** minute.  
**Note:** Store RNA at **-20°C**

