



Azma Elixir Pajoooh

Cat. No. AFTDX-1055

For Research Use Only

Fast **DNrich** Saliva Kit

Introduction

The Fast **DNrich** Saliva Kit provides all of the reagents necessary to extract DNA from saliva. DNA purified with this kit is suitable for a variety of applications, including amplification and digestion with restriction endonucleases.

Fast **DNrich** Saliva Kit components*

Cat. No. AFSDX1055-50

ST Solution	50 ml
TD Buffer	20 ml
Activator Reagent**	2 ml
VI Buffer	10 ml
SE Buffer	15 ml
Wash Buffer (conc.) **	50 ml
Elution Buffer	3 ml
Column	50
Manual	1

* 1.5 ml microcentrifuge tube, molecular biology grade water, and absolute ethanol are needed but are not included.

** Please refer to **Activator Reagent Preparation** and **Wash Buffer Preparation** before using this kit.

Chemical Hazard

Always wear gloves and practice standard safety precautions while using the kit. Do NOT disinfect extraction waste in solutions containing **bleach** or any other form of acid. To clean any items contaminated with the reagent, simply soak in detergent and water to remove all traces of contamination before cleaning with bleach or acidic solutions.

Reminder

Pre-set heater block at 65°C.

Prepare **Activator Reagent** immediately prior to use. Prepared Activator Reagent must be kept at -20°C.

Activator Reagent Preparation

Add 2 ml of molecular biology grade water to **Activator Reagent**, and vortex it well.

Note: Mark the check box on the bottle and write the date.

Note: For the best results, the prepared Activator Reagent should be used immediately. Prepared Activator Reagent can be stored for up to 3 months or 6 months at 4°C and -20°C, respectively.

Wash Buffer Preparation

Add 35 ml molecular biology grade **Absolute Ethanol** to **Wash Buffer** bottle before first use and mark the check box on it.



PROTOCOL

Step a: Sample Preparation

- a1 Add **2 ml** of saliva to a sterile 2 ml microcentrifuge tube then Centrifuge at **10000 g** for **5 minutes**, **Decant** the supernatant into waste.
- a2 Add **1000 µl** of **ST solution** to the sample tube and vortex it well.
- a3 Transfer **300 µl** of sample to a new 1.5 ml microcentrifuge tube.

Step b: Lysis

- b1 Add **400 µl** of **TD Buffer** and **40 µl** of **Activator Reagent** to the sample tube and vortex vigorously.
- b2 Incubate at **65°C** for **30 minutes**.
Note: During incubation time, vortex the sample tube every 10 minutes.
- b3 Incubate at **85°C** for **10 minutes** and vortex every 5 minutes.
- b4 **(Optional)** Add **200 µl** of **VI Buffer** to sample tube and vortex vigorously and then keep at RT for **5 minutes**.
- b5 Centrifuge at **11500 g** for **5 minutes**.
Note: Do not disturb the phases.
- b6 Carefully transfer about **300µl of supernatant** to a new 1.5 ml tube.
- b7 Add **300µl of SE Buffer** to the tube, invert for 5 times and keep at room temperature for **3 minutes** and then transfer all the sample to a **spin column**.
- b8 Centrifuge at **2000 g** for **2 minutes** and **discard** the flow through.

Step c: Washing

- c1 Add **500 µl** of **Wash Buffer** to the column, centrifuge at **8000 g** for **1 minutes** and **discard** the flow through.
- c2 **Repeat step c1.**

Step d: Column Drying

- d1 Centrifuge at **8000 g** for **1 minute**.
- d2 **Discard** the flow through and place the column into a new **1.5 ml** microcentrifuge tube.

Step e: DNA Elution

- e1 Add **50 µl** of **Elution Buffer** to the center of column and let stay at **RT** for **3 minutes**.
- e2 Centrifuge at **10000 g** for **2 minutes**.

Note: Store DNA at -20°C.

Short Protocol

Step a: Sample Preparation

- a1 Add **2 ml** of saliva to a sterile 2 ml microcentrifuge tube then Centrifuge at **10000 g** for **5 minutes**, **Decant** the supernatant into waste.
- a2 Add **1000 µl** of **ST solution** to the sample tube and vortex it well.
- a3 Transfer **300 µl** of sample to a new 1.5 ml microcentrifuge tube.

Step b: Lysis

- b1 Add **400 µl** of **TD Buffer** and **40 µl** of **Activator Reagent** to the sample tube and vortex vigorously.
- b2 Incubate at **65°C** for **30 minutes**.
- b3 Incubate at **85°C** for **10 minutes**.
- b4 Centrifuge at **11500 g** for **5 minutes**.
- b5 Carefully transfer about **300µl of supernatant** to a **spin column**.
- b6 Centrifuge at **2000 g** for **2 minutes** and **discard** the flow through.

Step c: Washing

- c1 Add **500 µl** of **Wash Buffer** to the column, centrifuge at **8000 g** for **1 minutes** and **discard** the flow through.
- c2 **Repeat step c1.**

Step d: Column Drying

- d1 Centrifuge at **8000 g** for **1 minute**.
- d2 **Discard** the flow through and place the column into a new **1.5 ml** microcentrifuge tube.

Step e: DNA Elution

- e1 Add **50 µl** of **Elution Buffer** to the center of column and let stay at **RT** for **3 minutes**.
- e2 Centrifuge at **10000 g** for **2 minutes**.





Azma Elixir Pajoooh

Cat. No. AFSDX1055

For Research Use Only

Fast **DNrich** Saliva Kit

لاندیپاز

