

Technical Data Sheet

Cetrimide Agar

Pseudomonas Selective Agar Base

acc. harm. EP/USP/JP

Ordering number: 1.05284.0500 / 1.05284.5000

A modification of the medium proposed by Brown and Lowbury (1965) for the isolation and differentiation of *Pseudomonas aeruginosa* from various materials.

This medium complies with the specifications given by the harmonized methods of EP, USP, JP for Microbial Examination of Non-sterile Products: Tests for Specified Microorganisms.

Mode of Action

The use of cetrimide (cetyltrimethylammonium bromide) was recommended by Lowbury (1951) and other authors; this compound largely inhibits the growth of the accompanying microbial flora. According to Lowbury and Collins (1955), a concentration of 0.3 g/l inhibits the accompanying organisms satisfactorily and minimizes interference with the growth of *Ps. aeruginosa*. The pigment production of *Pseudomonas aeruginosa* is not inhibited when grown on this medium.

Goto and Enomoto (1970) recommend the addition of 15 µg nalidixic acid/ml to improve the inhibition of the accompanying microbial flora.

Typical Composition

Peptone from Gelatin	20 g/l
MgCl ₂	1.4 g/l
K ₂ SO ₄	10 g/l
N-Cetyl-N,N,N-trimethylammoniumbromide (Cetrimide)	0.3 g/l
Agar-Agar	13.6 g/l
Glycerol	10 ml/l

Preparation

Suspend 45.3 g/l. Add 10 ml glycerol/l. Autoclave (15 min at 121 °C). Pour plates.

The appearance of the plates is turbid and light yellowish.

The pH value at 25 °C is in the range of 7.0-7.4.

Experimental Procedure and Evaluation

Inoculate by spreading the sample on the surface of the plates.

Incubation: *Pseudomonas aeruginosa* 18 h at 30-35 °C, others 72 h.

Pseudomonas aeruginosa colonies produce the bluish, non-fluorescent pyocyanin as well as the yellow-green pigment pyoverdin that fluoresces under UV light. Further tests should be performed for further identification (Hugh and Leifson 1953, Kovacs 1956, Thornley 1960, Bühlmann et al. 1961 etc).

Note: Beside *Pseudomonas aeruginosa* also *Pseudomonas putida* and *Pseudomonas fluorescens* are able to grow on Cetrimide Agar at 30-35 °C. The most important *Pseudomonads* can be pre-differentiated following the characteristics in the table below.

Phenotypic differentiation of the most important *Pseudomonads* ¹⁾

Characteristic	<i>P. aeruginosa</i>	<i>P. fluorescens</i>	<i>P. putida</i>	<i>P. stutzeri</i>
Oxidase	+	+	+	+
Growth at 41 °C	+	-	-	+/-
Pyoverdin („Fluorescein“)	+	+	+	-
Pyocyanin	+	-	-	-
Gelatinase	+	+	-	-

¹⁾ Derived from Bergey's Manual of Determinative Bacteriology (1994) Ninth Edition. Williams & Wilkins (+ = 90 % or more of the strains are positive; - = 90 % or more of the strains are negative)

Storage

The product can be used for sampling until the expiry date if stored upright, protected from light and properly sealed at +15 °C to +25 °C.

After first opening of the bottle the content can be used up to the expiry date when stored dry and tightly closed at +15 °C to +25 °C.

Disposal

Please mind the respective regulations for the disposal of used culture medium (e.g. autoclave for 20 min at 121 °C, disinfect, incinerate etc.).



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

EMD Millipore and the M logo are registered trademarks of Merck KGaA, Darmstadt, Germany. Lit. No. TN1351ENUS

Quality Control

Control Strains	ATCC #	Inoculum CFU	Incubation	Expected Results
<i>Pseudomonas aeruginosa</i>	9027	10 - 100	18 h at 30-35 °C	Recovery ≥ 50 %, yellow-green pigments
<i>Pseudomonas aeruginosa</i>	25668	10 - 100	18 h at 30-35 °C	Recovery ≥ 50 %, yellow-green pigments
<i>Pseudomonas aeruginosa</i>	27853	10 - 100	18 h at 30-35 °C	Recovery ≥ 50 %, yellow-green pigments
<i>Escherichia coli</i>	8739	> 10 ⁴	72 h at 30-35 °C	No growth
<i>Proteus mirabilis</i>	29906	> 10 ⁴	72 h at 30-35 °C	No growth
<i>Salmonella</i> Typhimurium	14028	> 10 ⁴	72 h at 30-35 °C	No growth
<i>Staphylococcus aureus</i>	6538	> 10 ⁴	72 h at 30-35 °C	No growth

Please refer to the actual batch related Certificate of Analysis.



Pseudomonas aeruginosa ATCC 9027



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

EMD Millipore and the M logo are registered trademarks of Merck KGaA, Darmstadt, Germany. Lit. No. TN1351ENUS

Literature

Brown, V.I. and Lowbury E.J.L. (1965): Use of improved cetrimide agar medium and other culture methods for *Pseudomonas aeruginosa*. J. Clin. Pathol. **18**: 752-756.

Buhlmann, X., Fischer, W.A. and Bruhn, J. (1961): Identification of a pyocyanogenic strains of *Pseudomonas aeruginosa*. J. Bact. **82**: 787-788.

DIN Deutsches Institut für Normung e.V. (1982): DIN 38411-8. Deutsche Einheitsverfahren zur Wasser-, Abwasser- und Schlammuntersuchung. Mikrobiologische Verfahren. Nachweis von *Pseudomonas aeruginosa*.

European Directorate for the Quality of Medicines and Healthcare. (2014): The European Pharmacopoeia. 8th Ed. Chapter 2.6.13 Microbiological examination of non-sterile products: Test for specified products. Strasbourg, France.

Goto, S. and Enomoto, S. (1970): Nalidixic acid cetrimide agar. A new selective plating medium for the selective isolation of *Pseudomonas aeruginosa*. Japan J. Microbiol. **14**: 65-72.

Hugh, R. and Leifson, E. (1953): The taxonomic significance of fermentative versus oxidative metabolism of carbohydrates by various gramnegative bacteria. J. Bact. **66**: 24-26.

Japanese Ministry of Health, Labour and Welfare. (2011): The Japanese Pharmacopoeia. 16th Ed. Chapter 4.05 Microbial Limit Test II. Microbiological examination of non-sterile products: Test for specified products. Japanese Ministry of Health, Labour and Welfare. Tokyo, Japan.

Lowbury, E.J.L. (1951): Improved culture methods for the detection of *Ps. pyocyanea*. J. Clin. Pathol. **4**: 66-72.

Lowbury, E.J.L. and Collins, A.G. (1955): The use of a new cetrimide product in a selective medium for *Pseudomonas pyocyanea*. J. Clin. Pathol. **8**: 47-48.

Thornley, M.J. (1960). The differentiation of *Pseudomonas* from other gram-negative bacteria on the basis of arginine metabolism. J. Appl. Bact. **23**: 37-52.

United States Pharmacopoeia 38 NF 33 (2015): <62> Microbiological examination of non-sterile products: Tests for specified microorganisms.

Ordering Information

Product	Cat. No.	Pack size
Cetrimide Agar <i>Pseudomonas</i> Selective Agar Base	1.05284.0500	500 g
Glycerol (about 87 %)	1.04094.0500	500 ml
UV Lamp (366 nm)	1.13203.0001	1 piece

EMD Millipore Corporation
290 Concord Road
Billerica, MA 01821, U.S.A.
mibio@emdmillipore.com
www.emdmillipore.com/biomonitoring

Find contact information for your country at:
www.emdmillipore.com/offices
For Technical Service, please visit:
www.emdmillipore.com/techservice



We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

EMD Millipore and the M logo are registered trademarks of Merck KGaA, Darmstadt, Germany. Lit. No. TN1351ENUS